

2 GAS PROFILE®

Created by MODULAR SERVICES COMPANY



The Profile® is a vertical headwall that makes the healing environment clinically functional for more complete and dependable care. With a unique gas pod design that houses two medical gases at one condensed level, the 2-gas Profile® provides a compact solution for bedside patient treatment.



MODULAR
SERVICES COMPANY

www.modularservices.com

Clinical Benefits

- Provides bedside positioning to limit patient interference
- 2 vertical tracks for equipment and storage management
- Up to 2 gases positioned at one ergonomic level
- Provides a clear path to essential services and equipment on the headwall
- Choose from a wide variety of accessories to support clinical needs



Product Features

- ① Vertical headwall positioned at bedside, pre-piped and wired with mounting studs attached (ships to job-site as one complete unit)
- ② Sleek Design: the Profile® is recessed into the wall, only protruding 1.25 inches from the wall
- ③ Available Services
 - 2 gases at one ergonomic level
 - Option for a third gas
 - Double Duplex Emergency Power Receptacle
 - Double Duplex Normal Power Receptacle
 - Option for 1 data port and 1 multi-pin connector
 - Provisions for nurse call and code blue button
- ④ Choose from a wide variety of high-pressure laminate for fascia
- ⑤ The vertical track, trim and gas pod available in Surf or Platinum Grey
- ⑥ Simple Installation: Accommodates a wall finished thickness of 1.25 inches