


PROPOSAL ONLY NOT FOR CONSTRUCTION

NOTES:

1. SERVICES TO BE REVERSED ON OPPOSITE HAND UNITS.
2. ALL UNITS ARE TO BE BUILT FOR MOUNTING INTO A 3-5/8" (SINGLE SIDED) OR 6" (DOUBLE SIDED) METAL STUD WALLS WITH 5/8" DRYWALL ON EACH SIDE.
3. PICTURE FRAME IS SHOWN IN UP POSITION. FRAME SLIDES DOWN AND OVER TO CONCEAL GAS SERVICES. FRAMED (34"W X 24"H) ARTWORK BY OTHERS MOUNTING KIT INCLUDED.
4. ARTWORKS CAN TO BE POWDER COATED TO MATCH FORMICA #923 "SURF".



LEGEND:

ACCESSORIES NOT SHOWN (PER UNIT)

REQ	SYM	DESCRIPTION	REQ	SYM	DESCRIPTION	REQ	DESCRIPTION
	O	OXYGEN OUTLET					
	A	AIR OUTLET					
	V	VACUUM OUTLET					



500 EAST BRITTON ROAD
OKLAHOMA CITY, OK 73114
www.modularservices.com

TOTAL QTY:
X

ARTWORKS® RECESSED-MOUNTED, 3 GAS
MODEL 6813-R & 6823-R

AS SHOWN:
X

DRAWN BY:
XXX

DRAWING: **4XXXX**

OPPOSITE:
X

DATE:
XX/XX/XX

PROJECT #:
XXXX

SALES ORDER #:
XXXXX

HOSPITAL
PROJECT
CITY, STATE

PROJECT