


PROPOSAL ONLY NOT FOR CONSTRUCTION

NOTES:

1. SERVICES TO BE REVERSED ON OPPOSITE HAND UNITS.
2. ALL UNITS ARE TO BE BUILT FOR MOUNTING INTO 3-5/8" (SINGLE SIDED) OR 6" (DOUBLE SIDED) METAL STUD WALLS WITH 5/8" DRYWALL.
3. PICTURE FRAME IS SHOWN IN UP POSITION. FRAME SLIDES DOWN AND OVER TO CONCEAL GAS SERVICES. FRAMED (45"W X 24"H) ARTWORK BY OTHERS MOUNTING KIT INCLUDED.
4. ARTWORKS CAN TO BE POWDER COATED TO MATCH FORMICA #923 "SURF".



LEGEND:

REQ	SYM	DESCRIPTION	REQ	SYM	DESCRIPTION
2	O	OXYGEN OUTLET			
1	A	AIR OUTLET			
2	V	VACUUM OUTLET			
2	SL	VACUUM BOTTLE SLIDE			
2	OS	OXYGEN SWIVEL ASSEMBLY			
1	AS	AIR SWIVEL ASSEMBLY			
2	VS	VACUUM SWIVEL ASSEMBLY			
2	FSL	FLIP DOWN VACUUM BOTTLE SLIDE			

ACCESSORIES NOT SHOWN (PER UNIT)

REQ	DESCRIPTION


 500 EAST BRITTON ROAD
 OKLAHOMA CITY, OK 73114
 www.modularservices.com

HOSPITAL
 PROJECT
 CITY, STATE

TOTAL QTY: **X**
 ARTWORKS® RECESSED-MOUNTED, 5 GAS
 MODEL 6815-R & 6825-R

AS SHOWN: X
 DRAWN BY: XXX
 DRAWING: **4XXXX**

OPPOSITE: X
 DATE: XX/XX/XX
 PROJECT #: XXXX
 SALES ORDER #: XXXXX