


PROPOSAL ONLY NOT FOR CONSTRUCTION


NOTES:

1. KYDEX COVER TO MATCH FORMICA #923 "SURF" OR PLATINUM GRAY.
2. FACE PLATES TO BE POWDER COATED "SURF" OR PLATINUM GRAY.




PARTITION WALL

PARTITION WALL



WIRING DIAGRAM



LEGEND:

ACCESSORIES NOT SHOWN (PER UNIT)

REQ	SYM	DESCRIPTION	REQ	SYM	DESCRIPTION	REQ	DESCRIPTION
2	N2	DUPLEX RECEPTACLE, 8300 IVORY 20A, 125V, H.G., NORMAL POWER					
1	E1	SINGLE RECEPTACLE, 8310 IVORY, 20A, 125V, H.G., NORMAL POWER					
2	E2	DUPLEX RECEPTACLE, 8300 RED 20A, 125V, H.G., EMER. POWER					
1	LN	LED NIGHT LIGHT, PASS & SEYMOUR, #NTLFULLWCC6, 120V, LOUVERED FACE					
1	MP	MULTI-PIN BED RECEPT. PROVISION					
1	DP	DATA OUTLET PROVISION					


**MODULAR SERVICES COMPANY\***  
 500 EAST BRITTON ROAD  
 OKLAHOMA CITY, OK 73114  
 www.modularservices.com

HOSPITAL PROJECT  
 CITY, STATE

TOTAL QTY: **X**  
**BED DOCKER**  
 AS SHOWN: X DRAWN BY: XXX DRAWING: **4XXXX**  
 OPPOSITE: X DATE: XX/XX/XX PROJECT #: XXXX SALES ORDER #: XXXXX