


FLATWALL

4500-S SERIES



MODULAR
SERVICES COMPANY

The Flatwall 4500-S Series is a vertically-oriented, surface-mounted headwall with a sleek, clean, linear design. Medical gas and electrical services are placed at convenient, ergonomic levels for ease of reach, without disturbing the patient.



FEATURES

- Surface-Mounted
- NeoTrack® angled equipment management tracks position accessories away from the headwall, providing a clear path to all services and mounted technology
- Optional forward-facing vertical tracks
- 3D Laminated steel panels
- Wide Variety of 3D laminate options
- Surf or Platinum Grey extrusion
- Compact device spacing
- Available in varying widths to accommodate acuity level needs
- Prefabricated to be piped and wired at finished ceiling
- Removable upper panel for quick access to all services
- Easy to clean with no exposed fasteners
- Easy to install
- Ships Fully Assembled



WWW.MODULARSERVICES.COM

ALL RIGHTS RESERVED
© 2020 Modular Services Company
DNC 70-1033, Rev. 7/20