


PROPOSAL ONLY NOT FOR CONSTRUCTION

NOTES:

1. SERVICES TO BE REVERSED ON OPPOSITE HAND UNITS.
2. ALL UNITS ARE TO BE BUILT FOR MOUNTING INTO 9" METAL STUD WALLS WITH 5/8" DRYWALL.
3. PICTURE FRAME IS SHOWN IN UP POSITION. FRAME SLIDES DOWN AND OVER TO CONCEAL GAS SERVICES. FRAMED (34"W X 30"H) ARTWORK BY OTHERS MOUNTING KIT INCLUDED.
4. ARTWORKS CAN TO BE POWDER COATED TO MATCH FORMICA #923 "SURF".



LEGEND:

ACCESSORIES NOT SHOWN (PER UNIT)

REQ	SYM	DESCRIPTION	REQ	SYM	DESCRIPTION	REQ	DESCRIPTION
1	O	OXYGEN OUTLET					
1	A	AIR OUTLET					
1	V	VACUUM OUTLET					
1	SL	VACUUM BOTTLE SLIDE					
2	N2	DUPLEX RECEPTACLE, 8300 IVORY 20A, 125V, H.G., NORMAL POWER					
2	E2	DUPLEX RECEPTACLE, 8300 RED 20A, 125V, H.G., EMER. POWER					



500 EAST BRITTON ROAD
OKLAHOMA CITY, OK 73114
www.modularservices.com

TOTAL QTY:
X

LARGE OPENING ARTWORKS®
MODEL 6830

AS SHOWN:
X

DRAWN BY:
XXX

DRAWING:

4XXXX

OPPOSITE:
X

DATE:
XX/XX/XX

PROJECT #:

XXXX

SALES ORDER #:

XXXXX

HOSPITAL
PROJECT
CITY, STATE

PROJECT