



LEGEND

REQ	SYM	DESCRIPTION
1	□	OXYGEN OUTLET
1	A	AIR OUTLET
1	V	VACUUM OUTLET
1	SL	VACUUM BOTTLE SLIDE
1	E2	DUPLEX RECEPTACLE, 8300 IVORY, 20A, 125V, H.G., NORMAL POWER
1	E4	DUPLEX RECEPTACLE, 8300 RED, 20A, 125V, H.G., EMER. POWER

QTY	ACCESSORIES NOT SHOWN (PER UNIT)

- NOTES:**
- ① FASCIA PANEL POWDER COATED TO MATCH 'SURF' OR 'PLATINUM GREY'. PERIPHERAL FRAMING IS CLEAR ANODIZED ALUMINUM.
 - ② GAS OUTLETS TO BE MSC STANDARD BRAND OUTLET WITH CONNECTION TO MATCH EXISTING HOSPITAL EQUIPMENT



**MODULAR SERVICES COMPANY**
500 EAST BRITTON ROAD
OKLAHOMA CITY, OK 73114


PROJECT

REVISED

DATE

DRAWN BY

FLUSH MOUNT CONSOLE

SCALE:  = 1'

DRAWING #
PB-CONS

— REQ'D: — AS SHOWN
— OPPOSITE

SERVICES TO BE REVERSED ON OPPOSITE HAND UNITS