

STRATUS[®]

Optimum Clinical Flexibility for Critical Care Environments



STRATUS[®]

The Stratus[®] was created to provide optimum flexibility for the positioning and accessibility of essential medical gases, electrical services and devices in critical care environments. The compact design gives the nursing staff complete access and clinical flexibility to quickly meet the needs of patients. Equipped with an Electromagnetic Braking System, the Stratus[®] service modules move quietly, smoothly, and position with ease. Stratus[®] is a dynamic solution for a practical need.



Service Module Options

5-Sided Column

10 gases at one ergonomic level
5 equipment management tracks
400 lbs. accessory load capacity

9x9 Column

8 gases at one ergonomic level
12 equipment management tracks
400 lbs. accessory load capacity

Stratus® Features

- 1 Arms have a 72-inch reach for complete flexibility to move the services and equipment around the bed and provide open access to the patient
- 2 Each articulating arm is capable of rotating up to 350 degrees at each joint for clinically correct positioning
- 3 Ergonomically placed medical gases and electrical services
- 4 Moving the Stratus® is smooth and fluid by using 3 Electromagnetic Brakes, controlled by an ergonomically placed handle
- 5 Vertical Tracks provide a place to mount essential equipment, accessories, shelving, and storage modules with supply drawers
- 6 Lighting Options: AIM LED Procedure Light Option and Ambient Lighting Option



Clinical Benefits

- Convenient movement to allow clinically correct positioning of essential services
- Open access to the patient's head
- Nurses can easily and quickly position all equipment, including monitors and pumps
- Ergonomically placed handle for easy movement
- Up to 10 medical gases at one level (5-sided Service Module)
- Vertical Tracks for equipment storage and supply management

CUSTOMIZED HEALING

Maximize ICU and ED patient care capabilities by combining other Modular products with the Stratus®. Choose from our extensive line of headwalls and columns to customize your space for an effective and efficient healing environment.

WHAT WE DO

Modular Services Company designs and manufactures advanced prefabricated products for the patient space: Medical Headwalls, Footwalls, Sinkwalls, Fixed/Rotating Columns, Movable Boom Systems and more. Our products are designed to be user-centric and meaningful. We strive to increase functionality for the clinical staff in the bedside environment. We offer smart, high design/aesthetic options for the design team and manufacture systems that make construction of the patient space both simple and efficient.

For more information about Modular Services Company and the products we create, please visit www.modularservices.com.

